



Spectral Wavelength Calibration Sources

Model Numbers: AR-2, HG-2, XE-2, KR-2, NE-2

NIST-traceable calibration of these wavelength calibration sources is not required. The wavelengths of the atomic emission lines produced by each source are determined by quantum mechanics and fixed as a function of time.

Therefore, these products do NOT require National Institute of Standards and Technology (NIST) certification.

The wavelength of the atomic emission lines may be verified at the NIST website:

http://physics.nist.gov/cgi-bin/AtData/lines_form

Nick Sebastian

Research & Science Sales Director, Americas

Ocean Optics, Inc.

p: 727-733-2447

a: 8060 Bryan Dairy Rd, Largo, FL 33777
4301 Metric Dr, Winter Park, FL 32792

w: www.oceanoptics.com e: nick.sebastian@oceanoptics.com



The Americas
830 Douglas Avenue
Dunedin, FL 34698
+1 727-733-2447
+1 727-733-3962 fax
www.oceanoptics.com
info@oceanoptics.com

Europe
Geograaf 24
6921 EW Dulven
The Netherlands
+31 26-319-0500
+31 26-319-0505 fax
www.oceanoptics.eu

Asia
666 Gubei Road, Kirin Tower
Suite 601B - Changning District
Shanghai, PRC 200336
+86 21-6295-6600
+86 21-6295-6708 fax
www.oceanopticschina.cn